

# VAHTS Universal Pro DNA Library Prep Kit for Illumina

Catalog# ND608



Version 9.1

Vazyme biotech co., ltd.

## 01/ Introduction

VAHTS Universal Pro DNA Library Prep Kit for Illumina is specially designed for library preparation for next generation sequencing (NGS) on Illumina® platforms, containing FFPE DNA repair modules, suitable for library preparation from 100 pg - 1 µg of input DNA. This kit contains a DNA damage repair module that can effectively repair DNA damage caused by formalin-fixed paraffin-embedded (FFPE), such as deamination of cytosine, nicks and gaps, oxidized bases, blocked 3' ends, etc, compatible with common DNA samples without affecting the quality of normal DNA sample libraries. Through the overall improvement of the module of end-repair, ligation and library amplification, the library conversion rate and the amplification output are greatly improved. It is widely applicable to PCR or PCR-Free library construction of multiple samples, and is compatible with targeted capture process.

## 02/ Components

Components	ND608-01 (24 rxn)	ND608-02 (96 rxn)
<input type="checkbox"/> DNA Damage Repair Enzyme	48 µl	192 µl
<input checked="" type="checkbox"/> End Prep Buffer	240 µl	960 µl
<input checked="" type="checkbox"/> End Prep Enzyme	120 µl	480 µl
<input checked="" type="checkbox"/> Rapid Ligation Buffer 2	600 µl	4 × 600 µl
<input checked="" type="checkbox"/> Rapid DNA Ligase	120 µl	480 µl
<input checked="" type="checkbox"/> VAHTS HiFi Amplification Mix	600 µl	4 × 600 µl
<input checked="" type="checkbox"/> PCR Primer Mix 3 for Illumina	120 µl	480 µl
<input checked="" type="checkbox"/> Control DNA (264 bp, 50 ng/µl)	10 µl	10 µl

## 03/ Storage

All components should be stored at -30°C ~ -15°C, and transported at -20°C ~ 0°C.

## 04/ Applications

Applicable to DNA library preparation for NGS on Illumina® platforms and compatible with various kinds of input samples, including Genomic DNA, cell-free DNA (cfDNA, ctDNA), formalin-fixed paraffin-embedded DNA (FFPE DNA), Chromatin immunoprecipitation DNA (ChIP DNA), and Amplicons. Input DNA is 100 pg - 1 µg. It is recommended to use this kit for:

- ◇ Whole genome sequencing.
- ◇ Whole exome or targeted sequencing (using Roche® NimbleGen™ SeqCap™ EZ, Agilent® SureSelect, or IDT xGen™ Lockdown™ Probes or other hybridization capture systems).
- ◇ Amplicon sequencing.
- ◇ ChIP-seq.
- ◇ Metagenome sequencing.
- ◇ Methylation Sequencing (in combination with Phanta UC Super-Fidelity DNA Polymerase for Library Amplification, Vazyme, #P507 and EpiArt DNA Methylation Bisulfite Kit, Vazyme #EM101)

## 05/ Additional Materials Required

- ◇ **Purification Beads:** VAHTS DNA Clean Beads (Vazyme #N411);
- ◇ **DNA quality control:** Agilent Technologies 2100 Bioanalyzer or other equivalent product;
- ◇ **DNA Adapters:**

VAHTS Dual Index UMI DNA Adapters for Illumina® (Vazyme #N331/332/333/334);

#N331-N334 are adapters with an 8 bp Unique Dual Index at two ends and a 10 bp UMI (Unique Molecular Identifier)

VAHTS Multiplex Oligoes Set 4/5 for Illumina® (Vazyme #N321/N322);

#N321/N322 are adapters with an 8 bp Unique Dual Index at two ends, allowing the construction of up to 384 different libraries by adapter combination.

VAHTS DNA Adapters set 1-6 for Illumina® (Vazyme #N801/N802 or #N805/N806/N807/N808);

#N801/N802 are 6 bp single-Index adapters, each containing 12 kinds. #N805/N806/N807/N808 are 8 bp-single-Index adapters, each containing 24 kinds.

- ◇ **Other Materials:**

Absolute ethanol, Ultrapure sterile water (ddH<sub>2</sub>O), 0.1 × TE, Elution Buffer (10 mM Tris-HCl, pH 8.0 - pH 8.5), Low adsorption EP tubes, PCR tubes, Magnetic stand, PCR instrument.



Vazyme Biotech Co., Ltd.  
www.vazyme.com

Order: global@vazyme.com

Support: global@vazyme.com

**For research use only, not for use in diagnostic procedures.**

## 06/ Notes

The parameters of library preparation procedures may be adjusted according to sample types, experimental designs, instruments, and operations. To obtain libraries of high quality, please read the following notes carefully.

For any questions during procedures, please contact Vazyme for help at [global@vazyme.com](mailto:global@vazyme.com).

### 06-1/ Input DNA & Fragmentation

◇ The recommended input DNA range is 100 pg - 1 µg. If possible, please use high quality DNA (with A260/A280 ratio of between 1.8 and 2.0) for library preparation. The recommended input DNA amounts are listed in **Table 1**.

**Table 1. Recommended Amount of Input DNA**

Application	Sample Type	Recommended Amount of Input DNA
Whole Genome Sequencing	Complex gDNA	50 ng - 1µg
Targeted Sequencing	Complex gDNA	10 ng - 1µg
Whole Genome/Targeted Sequencing	FFPE DNA	≥ 50 ng
Whole Genome/Targeted Sequencing	cfDNA/ctDNA	≥ 100 pg
Whole Genome Sequencing	Microbial genome	1 ng - 1µg
Whole Genome Sequencing (PCR-Free)	Complex/Simple genome	≥ 100 ng (no size selection) ≥ 200 ng (with size selection)
ChIP-seq	ChIP-seq	≥ 100 pg

▲ The amount recommended above is for DNA with high quality. For DNA with low quality, however, the input amount should be increased.

◇ "Input DNA" typically refers to the input DNA added to the [DNA Damage Repair & End Preparation](#), instead of the input amount for Fragmentation. If DNA was quantified before fragmentation, and fragmented DNA was subjected to cleanup or size selection prior to [DNA Damage Repair & End Preparation](#), the actual input into library construction may need to be recalculated. Otherwise, low library amplification cycle number may lead to low library yield.

◇ If proceeding with size selection, an elution volume of 105 µl is recommended. If proceeding directly to other subsequent steps, the recommended elution volume is 22.5 µl.

◇ DNA preparations containing high concentrations of EDTA, other chelating agents, or salts may affect the [DNA Damage Repair & End Preparation](#) reaction.

(1) Purify or size select fragmentation products and then dilute the purification products to 0.1× TE or ddH<sub>2</sub>O (≤ 48 µl) before library construction.

(2) If fragmented DNA is not subjected to a bead-based cleanup or size selection prior to library construction:

For mechanical fragmentation, DNA should be fragmented in 0.1× TE. Fragmentation in ultrapure sterile water is not recommended.

For enzymatic fragmentation, make sure there is not an excess of EDTA or other chelating agents in the [Stop Buffer](#).

### 06-2/ Adapters

◇ For the Illumina sequencing platform, Vazyme offers four sets of Indexed Adapters, which can be selected according to different usage needs and the number of Pooling samples:

▲ Vazyme, # N331/N332/N333/N334: up to 96 kinds of non-repeating dual-ended 8 bp Indexed adapters with 10 bp UMI sequences, 24 kinds/each set.

▲ Vazyme, #N321/N322: up to 384 kinds of dual-ended 8 bp Indexed Adapters.

▲ Vazyme, #N801/N802: up to 24 kinds of single-ended 6 bp Indexed Adapters, 12 kinds/each set;

▲ Vazyme, #N805/N806/N807/N808: up to 96 kinds of single-ended 8 bp Indexed Adapters, 24 kinds/each

◇ The quality and amount of Adapters directly affect the preparation efficiency and library quality. The recommended ratio of adapter: input is between 10 : 1 and 200 : 1. High Adapter input may lead to residual Adapter/Adapter Dimer. Low Adapter input may affect ligation efficiency and reduces library yields. Please refer to **Table 2** for the recommended adapter concentrations for different DNA inputs.

**Table 2. Recommended adapter concentrations for libraries prepared from 100 pg - 1 µg input DNA**

Input DNA	Adapter: Input DNA Molar Ratio	Adapter concentration from other source (Working concentration)	Vazyme Adapter Dilution Ratio
500 ng - 1 µg	10:1 - 20:1	10 µM	Undiluted
100 ng - 500 ng	20:1 - 100:1	10 µM	Undiluted
25 ng - 100 ng	50:1 - 200:1	5 µM	1:2
5 ng - 25 ng	40:1 - 200:1	1 µM	1:10
100 pg - 5 ng	60:1 - 3000:1	0.2 µM	1:30 - 1:200

▲ Calculate the moles of Input DNA:

moles of Input DNA (pmol) ≈ mass of Input DNA (ng) / [0.66 × average length of Input DNA (bp)]

▲ According to the concentration or dilution ratio, dilute Adapter with 0.1× TE. Make the volume of Adapter fixed (5 µl) to avoid pipetting error.



## 06/ Notes

▲ The quality of adapters will affect the molar ratio of Adapter and Input DNA and further affect ligation rate and library yields. Please see the adapter with high quality for library preparation. Dilute and store the Adapter solution with 0.1× TE. Minimize the number of freeze-thaw cycles.

▲ Increasing adapter inputs can increase library yields, especially when the amount of Input DNA is ≤ 25 ng. When optimizing workflows to increase the efficiency of library construction, several higher adapter concentrations should be evaluated: based on the recommended adapter concentration (Table 2), try several additional concentrations in a range that is 2-10 times higher than the recommended concentration. If the adapter concentration is limited, try using more volume to increase adapter amount. For example, if Input DNA is 500 ng-1 μg while default volume of adapter is 5 μl, please increase to 10 μl to enhance 5% - 15% library output. However, it should be noted that increasing the concentration of adapter may increase the residue of adapter in the library, resulting in wasted sequencing data.

### 06-3/ Cleanup of Adapter Ligation Products

◇ Unused Adapters should be removed before library amplification (for PCR amplification library) or sequencing (for PCR-free library). The default purification condition 0.6× (60 μl beads/100 μl products) is suitable for most cases. To obtain libraries with larger insert sizes, the amount of beads can be reduced to lower the content of small DNA fragments. Please note this is just a rough adjust. To control the library distribution accurately, please process size selection after cleanup.

◇ If proceeding with size selection, an elution volume of 105 μl is recommended. If proceeding directly to other subsequent steps, the recommended elution volume is 22.5 μl.

◇ A second cleanup may be performed using a 1× bead to DNA ratio, if post-ligation analysis reveals unacceptable levels of adapter and/or adapter-dimer carry-over after the first cleanup. Make the volume of the purification products from the first round up to 50 μl with ddH<sub>2</sub>O, and then add 50 μl of beads for second round purification. A second cleanup may be particularly beneficial when libraries are prepared in PCR-free workflows for direct sequencing on Illumina platforms. Sometimes it may be necessary to reduce the amount of adapter to eliminate the residual adapter and/or adapter-dimer completely.

### 06-4/ Beads

◇ This protocol has been validated for use with VAHTS DNA Clean Beads (Vazyme, #N411).

**Note:** If you use beads from other vendors, the purification conditions may need to be changed.

◇ General notes on beads manipulations:

▲ The amount of beads is calculated by "×" (multiple), indicating the multiple of beads volume compared to sample volume. For example, if sample volume is 100 μl, 1× beads means the volume of beads is 1 × 100 μl = 100 μl; 0.6×/0.2× size selection means the first round of bead volume is 0.6 × 100 μl = 60 μl and the second round is 0.2 × 100 μl = 20 μl.

▲ The volume of beads directly affects the purified DNA size of lower limit. The higher multiple, the smaller insert of lower limit, and vice-versa. For example, 1× beads can only purify DNA longer than 250 bp. The smaller fragments will be discarded during cleanup. While 1.8× beads can purify DNA of 150 bp.

▲ Equilibration to room temperature (place in room temperature for 30 min) before use is essential to achieve specified size distribution and high yield of libraries.

▲ Beads will settle gradually. Always ensure that they are fully resuspended before use by vortexing or up-and-down pipetting several times.

▲ The time required for complete capture of beads varies according to the reaction vessel and magnet used. It is important not to discard or transfer any beads with the removal or transfer of supernatant. Please transfer the supernatant after the solution is completely clarified, and leave 2 μl - 3 μl of supernatant behind to avoid disturbing the beads. If the magnetic beads are accidentally drawn out with the supernatant, the yield will decrease, the effect of size selection will be poor, and even the subsequent enzymatic reaction will be affected. In this case, the magnetic beads can be mixed and placed on magnetic stand again to separate the beads completely.

▲ Always use freshly prepared 80% ethanol. Keeping tubes on magnet stand without disturbing the beads during elution.

▲ It is important to remove all the ethanol before proceeding with subsequent reactions. However, over-drying of beads may make them difficult to resuspend, resulting in a dramatic loss of DNA. Normally, drying of beads for 5 min–10 min at room temperature should be sufficient. Do not heat and dry (such as oven drying at 37°C).

▲ DNA should be eluted from beads with elution buffer (10 mM Tris-HCl, pH 8.0 - pH 8.5) for stable preservation. However, DNA must be eluted and stored in PCR-grade water (ultrapure sterile water) for subsequent target capture reaction, to facilitate concentrating of DNA libraries before probe hybridization.

▲ Purified DNA in elution buffer should be stable at 4°C for 1 week, or at -20°C for long-term storage. Avoid excessive freezing and thawing cycles.

### 06-5/ Size Selection

◇ If the distribution range of input DNA is broad, size selection will be necessary to control the final library size. It is recommended to use a dual bead-based size selection, while gel-based size selection technique is also usable.

◇ Size selection may be carried out at several time points in the overall workflow, for example:

▲ After dsDNA Fragmentation;

▲ after the cleanup of adapter ligation products;

▲ after library amplification.



## 06/ Notes

◇The standard protocol of this manual (Refer to **08/ Standard Protocol for Library Preparation**) does not include size-selection. Please refer to **Appendix 1** for detailed protocols for size-selection.

◇Size selection inevitably leads to a loss of sample material. These losses can be 60% - 95%. The potential advantages of size selection in a library construction workflow should be weighed against the potential loss of library complexity, especially when input DNA is limited. Two or more size selection steps will result in dramatic decrease in library complexity and yields.

◇Over-amplification typically results in the observation of secondary, higher molecular weight peaks in the electrophoretic profiles of amplified libraries. These higher molecular weight peaks are artifacts (Refer to **06-6/ Library Amplification**) of the analysis, and typically contain authentic library molecules of the appropriate length. To eliminate these artifacts, optimization of library amplification reaction parameters (cycle number and primer concentration), rather than post-amplification size selection, is recommended.

◇**Rapid Ligation Buffer 2** contains a high concentration of PEG, which, if not removed, will interfere with efficient dual bead-based size selection and gel-based size selection. Therefore, if size selection is performed after Adapter Ligation, it is important to perform at least one step of bead-based clean-up (Refer to **08/ Standard Protocol for Library Preparation, Step 02. Adapter Ligation, Step 6. Purify the Adapter Ligation products using VAHTS DNA Clean Beads**) prior to performing bead- or electrophoresis-based size selection. If size selection is performed after Library Amplification, the original clean-up step can be directly replaced by dual bead-based size selection or gel-based size selection.

### 06-6/ Library Amplification

◇The PCR Primer Mix 3 is suitable for the amplification of all Illumina® libraries flanked by the P5 and P7 flow cell sequences. User-supplied primer mixes may be used in combination with incomplete or custom adapters. Each primer should be used at a final concentration of 5 µM-20 µM.

◇In library amplification reactions, primers are typically depleted before dNTPs. When DNA synthesis can no longer take place due to substrate depletion, subsequent rounds of DNA denaturation and annealing result in the formation of improperly annealed, partially double-stranded, heteroduplex DNA. These species migrate slower and are observed as secondary, higher molecular weight peaks during the electrophoretic analysis of amplified libraries. However, they typically comprise library molecules of the desired length, which are individualized during denaturation prior to cluster amplification or probe hybridization. Since these heteroduplexes contain significant portions of single-stranded DNA, over-amplification leads to the under-quantification of library molecules with assays employing dsDNA-binding dyes (i.e. **Euqalbit dsDNA HS Assay Kit, Vazyme, #EQ111**). qPCR-based library quantification methods, such as the **VAHTS Library Quantification Kit for Illumina® (Vazyme, #NQ101-NQ106)**, quantify DNA by denaturation and amplification, thereby providing an accurate measure of the amount of adapter-ligated molecules in a library - even if the library was over-amplified.

◇The number of amplification cycles for library amplification should be limited as much as possible. Insufficient library amplification leads to insufficient library output. Excessive library amplification can result in other unwanted artifacts such as amplification bias, PCR duplicates, chimeric library inserts and amplified mutations. **Table 3** provides recommended cycle numbers for libraries prepared from 100 pg-1 µg high-quality input DNA, to obtain approximately 100 ng or 1 µg of amplified library.

**Table 3. Recommended cycle numbers for 100 pg–1 µg of input DNA**

Input DNA (Into End Preparation)	Number of cycles required to generate	
	100 ng	1 µg
100 pg	13 - 16	16 - 19
1 ng	9 - 11	12 - 15
5 ng	6 - 9	9 - 13
10 ng	5 - 8	8 - 12
50 ng	3 - 4	6 - 9
100 ng	2 - 3	5 - 8
250 ng	1 - 2	3 - 6
500 ng	0	2 - 5
1 µg	0	1 - 3

▲ The table above is the recommended cycle numbers for high-quality input DNA of 200 bp. The quality of FFPE DNA varies greatly. If the input DNA quality is low or library length is long, please increase cycle numbers to obtain a sufficient library.

▲ If size selection is proceeded after adapter ligation, please choose larger cycle numbers for Library Amplification. Otherwise, please choose the small numbers.

◇If the adapters are complete (e.g. **Vazyme, #N801/N802/N805/ N806/ N807/ N808 or #N331/N332/N333/N334**) and a sufficient amount of library is available, it may be possible to skip library amplification to obtain PCR-Free libraries.

◇When using incomplete adapters (e.g. **Vazyme, #N321/N322**), a minimum number of amplification cycles (at least 2) is required to complete adapter sequences


## 06/ Notes

### 06-7/ Evaluation of Library Quality

Normally, a constructed library can be evaluated for quality by size distribution and concentration detection.

#### ◇Size Distribution

▲The size distribution of final libraries can be confirmed with an electrophoretic-based method. A LabChip® GX, GXII or GX Touch (PerkinElmer), Bioanalyzer or TapeStation (Agilent Technologies), Fragment Analyzer™ (Advanced Analytical) or similar instrument is recommended. Typical electrophoretic profiles for libraries prepared with VAHTS Universal Pro DNA Library Prep Kit for Illumina are shown in Fig. 1.



**Fig. 1. Length Distribution of Products during Library Preparation using VAHTS Universal Pro DNA Library Prep Kit for Illumina**

▲Please note that libraries prepared with “forked” adapters in PCR-free workflows will appear to have a longer than expected mode fragment length, and/or may display a broad or bimodal size distribution when analyzed electrophoretically. To accurately determine the size distribution of an unamplified library, an aliquot of the library may be subjected to a few cycles of amplification prior to electrophoretic analysis, to ensure that all adapter-ligated molecules are fully double-stranded. Alternatively, size information may be obtained by electrophoretic analysis of library quantification products generated with VAHTS Library Quantification Kit (Vazyme, #NQ101-NQ106).

#### ◇Quantification of Libraries

There are two methods of library quantification: one is based on dsDNA fluorescent dyes, i.e. Qubit®, PicoGreen®, or Equalbit dsDNA HS Assay Kit (Vazyme, #EQ111); the other is qPCR-based quantification, such as VAHTS Library Quantification Kit (Vazyme, #NQ101-NQ106). The first one is easy to proceed, however, qPCR-based method is recommended due to the following reasons:

▲When using full-length adapters, and once ligation has been completed, qPCR-based quantification kit can quantify libraries at different stages of the workflow. Thus, the efficiency of End Preparation, purification/size selection, and Library Amplification can be assessed, to provide useful data for optimization or troubleshooting.

▲PCR-Free libraries contain some fragments with single-end adapters or without adapters. When using the method of double-stranded DNA dye (Equalbit dsDNA HS Assay Kit, Vazyme, #EQ111), these fragments will be also measured. But qPCR quantification only quantifies those molecules with double-end adapters in the correct orientation for sequencing. Therefore, PCR-Free libraries can only be quantified by qPCR-based method.

▲Over-amplified libraries contain non-complete fragment and can't be measured by Qubit® Equalbit dsDNA HS Assay Kit (Vazyme, #EQ111) or PicoGreen® methods. Measurements with qPCR-based method are not affected by library over amplification.

### 06-8/ Other Notes

◇Thaw all the components at room temperature before use. Mix thoroughly by turning up and down multiple times after thawing. Centrifuge briefly and place on ice.

◇Mixing with pipetting is recommended when preparing solutions. Shaking excessively may lead to reduced yield.


◇To avoid cross contamination, tips with filter are recommended. Change tips between samples.

◇It is recommended to use PCR instrument with heated lid. Preheat PCR instrument to reaction temperature in advance.

◇Aerosol contamination is easily to occur due to improper PCR operations, which affects experiment accuracy. Therefore, it is recommended to separate preparation area and clean-up area physically, use dedicated pipettor, and clean experimental region by 0.5% sodium hypochlorite or 10% decolorizer timely.

07/ Mechanism & Workflow


07-1/ Mechanism



Mechanism of VAHTS Universal Pro DNA Library Prep Kit for Illumina

## 07/ Mechanism & Workflow

### 07-2/ Workflow



Workflow of VAHTS Universal Pro DNA Library Prep Kit for Illumina

## 08/ Standard Protocol for Library Preparation

### Step 01. DNA Damage Repair & End Preparation

This step is for DNA Damage Repair (caused by formalin-fixed paraffin-embedded (FFPE), such as nicks and gaps, oxidized bases), End Repair, 5' phosphorylated, and dA-tailing.

1. Thaw the End Prep Buffer and spin down briefly. Prepare the reaction solution in a PCR tube as follows:

Components	Volume	
Input DNA	x µl	
DNA Damage Repair Enzyme*	2 µl	<input type="checkbox"/>
End Prep Enzyme	5 µl	<input checked="" type="checkbox"/>
End Prep Buffer	10 µl	<input checked="" type="checkbox"/>
ddH <sub>2</sub> O	To 65 µl	

\* If the sample does not require DNA Damage Repair, do not add this enzyme.

2. Mix thoroughly by gently pipetting up and down. DO NOT Vortex! Spin down briefly.

## 08/ Standard Protocol for Library Preparation

3. Put the tube in a PCR instrument and run the following PCR program:

Temperature	Time
Hot lid of 105°C	On
30°C	30 min
65°C	15 min
4°C	Hold

### Step 02. Adapter Ligation

This step is to ligate the product of [DNA Damage Repair & End Preparation](#) from Step 01 to the adapter.

1. Dilute the Adapter to the appropriate concentration according to the amount of Input DNA, please refer to Table 2 (Page 04).

2. Thaw the Rapid Ligation buffer 2 and mix thoroughly. Place on ice.

3. Prepare the reaction solution in a sterile PCR tube as follows:

Components	Volume
End Preparation Products	65 µl
Rapid Ligation buffer 2	25 µl
Rapid DNA ligase	5 µl
DNA Adapter X	5 µl
Total	100 µl

▲ If VAHTS Multiplex Oligos Set 4/5 for Illumina (Vazyme #N321/N322) is used, the adapter should be of the DNA adapter-S for Illumina in the kit, and still be used in an amount of 5 µl.

4. Mix thoroughly by gently pipetting up and down. DO NOT Vortex! Spin down briefly.

5. Put the tube in a PCR instrument and run the following PCR program:

Temperature	Time
Hot lid of 105°C	On
20°C	15 min
4°C	Hold

▲ For the low-amount of Input DNA, try to double the ligation time. However, extending the ligation time may result in increase of the Adapter Dime. If necessary, the Adapter concentration should be adjusted at the same time.

6. Purify the Adapter Ligation products using VAHTS DNA Clean Beads:

1/ Equilibrate the VAHTS DNA Clean Beads to room temperature. Suspend the beads thoroughly by vortexing.

2/ Pipet 60 µl of beads into 100 µl of the Adapter Ligation products. Mix thoroughly by vortexing or pipetting up and down for 10 times.

3/ Incubate at room temperature for 5 min.

4/ Place the sample on a magnetic stand. Wait until the solution clarifies (about 5 min). Keep it on magnetic stand and carefully discard the supernatant without disturbing the beads.

5/ Keeping the sample on magnetic stand, add 200 µl of freshly prepared 80% ethanol to rinse the beads. **DO NOT re-suspend the beads!** Incubate for 30 sec at room temperature and carefully discard the supernatant without disturbing the beads.

6/ Repeat the **Step 5/**.

7/ Keep the tube on the magnetic stand, open the tube, and air-dry the beads for 5 - 10 min (avoid over-drying) until there is no ethanol residue.

8/ Take the tube out of the magnetic stand for elution:

▲ For products with no need for size-selection: Add 22.5 µl of elution buffer (10 mM Tris-HCl, pH 8.0 - pH 8.5) to elute DNA. Mix thoroughly by vortexing or pipetting and incubate for 2 min at room temperature. Place the tube back on the magnetic stand and wait until the solution clarifies (about 5 min). Carefully transfer 20 µl of the supernatant to a new Nuclease-free PCR tube without disturbing the beads.

▲ For products with need for size-selection: Add 105 µl of elution buffer (10 mM Tris-HCl, pH 8.0 - pH 8.5) to elute DNA. Mix thoroughly by vortexing or pipetting and incubate for 2 min at room temperature. Place the tube back on the magnetic stand and wait until the solution clarifies (about 5 min). Carefully transfer 100 µl of the supernatant to a new Nuclease-free PCR tube without disturbing the beads. Proceed to size-selection according to [Table 5](#) (Page 10).



The products can be stable for one week at 4°C. Keep at -20°C for long-term storage. Avoid unnecessary freeze-and-thaw cycles.

## 08/ Standard Protocol for Library Preparation

### Step 03. Library Amplification

This step is to amplify the purified or size-selected adapter ligation products. Whether to proceed with this step depends on the amount of input DNA, whether adapters are in complete length, and downstream application. If adapters are not in complete length (e.g., Vazyme, # N321/# N322), this step is necessary. If adapters are in complete length, for input DNA < 50 ng, library amplification is recommended. Skip this step and proceed directly to **Step 04**, if input DNA is ≥ 50 ng or there is no need for library amplification.

1. Thaw PCR Primer Mix 3 and VAHTS HiFi Amplification Mix, and mix thoroughly. Prepare the reaction solution in a sterile PCR tube as follows:

Components	Volume
Purified or size-selected adapter ligation products	20 µl
PCR Primer Mix 3 for Illumina	5 µl 
VAHTS HiFi Amplification Mix	25 µl 
Total	50 µl

▲ If VAHTS Multiplex Oligos Set 4/5 for Illumina (Vazyme #N321/N322) is used, the adapter should be of the DNA adapter-S for Illumina in the kit, and still be used in an amount of 5 µl.

2. Mix thoroughly by gently pipetting up and down. DO NOT Vortex! Spin down briefly.

3. Put the tube in a PCR instrument and run the following PCR program (Hot Lid Temperature: 105°C):

Temperature	Time	Cycles
95°C	3 min	1
98°C	20 sec	According to <b>Table 3</b> (Page 04)
60°C	15 sec	
72°C	30 sec	
72°C	5 min	1
4°C	Hold	

4. For size selection, please refer to **Appendix 1**. If there is no need for size selection, please purify the products with VAHTS DNA Clean Beads as follows:

1/ Equilibrate the VAHTS DNA Clean Beads to room temperature. Suspend the beads thoroughly by vortexing.

2/ Pipet 45 µl (0.9x) of beads into 50 µl of the Library Amplification products. Mix thoroughly by vortexing or pipetting up and down for 10 times.

3/ Incubate at room temperature for 5 min.

4/ Place the tube on a magnetic stand. Wait until the solution clarifies (about 5 min). Keep it on the magnetic stand and carefully discard the supernatant without disturbing the beads.

5/ Keeping the sample on the magnetic stand, add 200 µl of freshly prepared 80% ethanol to rinse the beads. **DO NOT re-suspend the beads!** Incubate for 30 sec at room temperature and carefully discard the supernatant without disturbing the beads.

6/ Repeat the **Step 5/**.

7/ Keep the tube on the magnetic stand, open the tube and air-dry the beads for 5 - 10 min.

8/ Take the Tube out of the magnetic stand for elution:

▲ For products with no need for Targeted Capture: Add 22.5 µl of elution buffer (10 mM Tris-HCl, pH 8.0 - pH 8.5) to elute DNA. Mix thoroughly by vortexing or pipetting and incubate for 2 min at room temperature. Place the tube back on the magnetic stand and wait until the solution clarifies (about 5 min). Carefully transfer 20 µl of the supernatant to a new Nuclease-free PCR tube without disturbing the beads.

▲ For products with need for Targeted Capture: Add 22.5 µl of ddH<sub>2</sub>O to elute DNA. Mix thoroughly by vortexing or pipetting and incubate for 2 min at room temperature. Place the tube back on the magnetic stand and wait until the solution clarifies (about 5 min). Carefully transfer 20 µl of the supernatant to a new Nuclease-free PCR tube without disturbing the beads.



The products can be stable for one week at 4°C. Keep at -20°C for long-term storage. Avoid unnecessary freeze-and-thaw cycles.

### Step 04. Quality Control of Library

Please refer to **06-7. Evaluation of Library Quality** for details.

## Appendix 1. Dual bead-based size selection

◇ To meet the needs of different sequencing applications, dual bead-based size-selection is necessary to control the distribution of Insert Size. Generally, the size-selection is recommended to be carried out after the post-ligation cleanup. It also can be arranged prior to End Preparation or after library amplification. Ensure the point of size selection is unique, for two or more times of size selection will lead to dramatic decrease in library complexity and yields. Please refer to Table 4 for the choice of size selection points and advantages/disadvantages of each point.

**Table 4. The Choice of Size Selection Point**

Size Selection Point	Applicable Conditions	Advantages	Disadvantages	Applicable Samples
After dsDNA Fragmentation	Sufficient input DNA, with broad size distribution or unexpected insert size; input DNA with low purity.	Selected length is concentrated; accurate amount of input DNA; increased DNA purity to increase library prep success rate;	Loss of DNA; size distribution is broad*	gDNA with insufficient or excess fragmentation
After Adapter Ligation	Sufficient input DNA with proper size distribution**	Decrease the loss of short input-DNA	Size distribution is broad*	Proper fragmentation of gDNA or FFPE DNA with broad size distribution
After Library Amplification	Low input DNA**	Decrease the loss of input DNA during workflow, increase library complexity	Size distribution is broad	cfDNA
No Size-selection	Size distribution of input DNA is proper; low input DNA	Decrease the loss of input DNA during workflow, increase library complexity	Can't control the insert size	Multiple PCR products; fragmented FFPE DNA

\* DNA ends affect the result of size selection. Single-stranded portion at the end of the Input DNA and single-stranded non-complementary arms of "Y" adapters ligated to DNA fragments may result in a broader size distribution.

\*\* If the amount of input DNA is  $\geq 100$  ng, size selection after Adapter Ligation is recommended. If the amount of input DNA is  $< 100$  ng, or the copies of samples is limited, size selection after Library Amplification is recommended.

◇ Dual bead-based size selection is used to select DNA fragments of expected length by controlling the amount of beads used. The mechanism of size-selection is: (1) during the first round, DNA with larger fragments bind to beads and are discarded with these beads (2) during the second round, DNA of expected length bind to bead while DNA with smaller fragments in the supernatant are discarded. Components in the initial DNA may affect the result of size selection. Therefore, according to the point of size selection, the amount of beads differs. Please refer to **Table 5** to choose the volume of beads according to expected insert size and selection points.

**Table 5. Size Selection of Library**

Point of size selection	Amount of Beads Added	Expected Insert Size (bp)									
		150	200	250	300	350	400	450	500	550	700
After dsDNA Fragmentation (Sample volume is added to 100 $\mu$ l)	1st-Round X( $\mu$ l)	100	90	80	70	60	55	52	50	48	43
	2nd-Round Y( $\mu$ l)	30	20	20	20	20	20	15	15	15	12
After Adapter Ligation (Sample volume is 100 $\mu$ l)	1st-Round X( $\mu$ l)	78	68	65	59	56	53	51	50	/	/
	2nd-Round Y( $\mu$ l)	20	20	15	15	12	12	10	10	/	/
After Library Amplification (Sample volume is adjusted to 100 $\mu$ l)	1st-Round X( $\mu$ l)	78	70	63	55	50	46	45	44	/	/
	2nd-Round Y( $\mu$ l)	20	20	20	20	20	20	20	15	/	/

◇ If adapters are not in complete length (e.g., Vazyme, # N321 / N322 ), please refer to the following table to choose the volume of beads according to expected insert size and selection points.

Point of size selection	Amount of Beads Added	Expected Insert Size(bp)							
		150	200	250	300	350	400	450	500
After Adapter Ligation (Sample volume is 100 $\mu$ l)	1st-Round X( $\mu$ l)	100	90	75	65	60	55	53	50
	2nd-Round Y( $\mu$ l)	20	20	20	20	20	20	20	18

▲ In bead-based size selection, the larger Insert Size, the broader size distribution. However, beads cannot select DNA with insert size  $> 700$  bp, please choose electrophoretic method instead.

▲ The volume ratio of samples and beads is important for size selection. Please ensure the accuracy of initial sample volume and pipetting volume.



Vazyme Biotech Co., Ltd.  
www.vazyme.com

Order: global@vazyme.com

Support: global@vazyme.com

**For research use only, not for use in diagnostic procedures.**

## Appendix 1. Dual bead-based size selection

### ◇ Pretreatment of Samples (Important!)

▲ For size-selection after dsDNA Fragmentation, Adapter Ligation purification or Library Amplification, the initial sample volume should be 100 µl. If the volume is less than 100 µl, please add to 100 µl using ddH<sub>2</sub>O.

▲ If no such pretreatment is performed, please adjust the beads volume according to the actual sample volume. However, small volume of samples will lead to increasing inaccuracy in pipetting, resulting in inaccuracy in size selection. Therefore, it is NOT recommended to size-select samples of < 50 µl directly without pretreatment.

### ◇ Protocol for Size Selection [Refer to **Table 5** (Page 10 ) for the value of X and Y]

1. Equilibrate the VAHTS DNA Clean Beads to room temperature. Suspend the beads thoroughly by vortexing.
2. Add X µl of beads into 100µl of the DNA products. Mix thoroughly by vortexing or pipetting up and down for 10 times.
3. Incubate at room temperature for 5 min.
4. Place the sample on the magnetic stand. Wait until the solution clarifies (about 5 min). Keep it on the magnetic stand and carefully transfer the supernatant to a new Nuclease-free PCR tube without disturbing the beads.
5. Add Y µl of beads into the supernatant from Step 4. Mix thoroughly by vortexing or pipetting up and down for 10 times.
6. Incubate at room temperature for 5 min.
7. Place the sample on the magnetic stand. Wait until the solution clarifies (about 5 min). Keep it on the magnetic separation rack and carefully discard the supernatant without disturbing the bead.
8. Keep the tube on the magnetic stand, add 200 µl of freshly prepared 80% ethanol to the beads. DO NOT re-suspend the beads! Incubate at room temperature for 30 sec and carefully discard the supernatant without disturbing the beads.
9. Repeat the Step 8.
10. Keep the tube on the magnetic stand, open the tube and air-dry the beads for 5 - 10 min.
11. Take the Tube out of the magnetic stand for elution.

▲ For products with no need for Targeted Capture: Add 22.5 µl of elution buffer (10 mM Tris-HCl, pH 8.0-pH 8.5) to elute DNA. Mix thoroughly by vortexing or pipetting and incubate for 2 min at room temperature. Place the tube back on the magnetic stand and wait until the solution clarifies (about 5 min). Carefully transfer 20 µl of the supernatant to a new Nuclease-free PCR tube without disturbing the beads

▲ For products with need for Targeted Capture: Add 22.5 µl of ddH<sub>2</sub>O to elute DNA. Mix thoroughly by vortexing or pipetting and incubate for 2 min at room temperature. Place the tube back on the magnetic stand and wait until the solution clarifies (about 5 min). Carefully transfer 20 µl of the supernatant to a new Nuclease-free PCR tube without disturbing the beads.

## Appendix 2. cfDNA Library Preparation

Cell-free DNA (cfDNA) is a highly fragmented (~180 bp) DNA in blood with low content. cfDNA is of great importance in noninvasive prenatal testing (NIPT) and liquid biopsy. VAHTS Universal Pro DNA Library Prep Kit for Illumina provides an easy and optimized solution for the library preparation of cfDNA.

### Notes

- ◇ "Input DNA" is the DNA added to DNA Damage Repair & End Preparation with volume ≤ 48 µl.
- ◇ cfDNA is highly fragmented. No fragmentation is needed for cfDNA.
- ◇ To ensure the quality of library, it is recommended to detect the size distribution (by an Agilent 2100 Bioanalyzer) and concentration (by Equalbit) of input DNA.

### Library Preparation Procedures

#### Step 1: DNA Damage Repair & End Preparation (Refer to 08/ Standard Protocol of Library Preparation, Step 01)

Input DNA: 100 pg-100 ng, this step does not use DNA Damage Repair Enzyme.

#### Step 2: Adapter Ligation (Refer to 08/ Standard Protocol of Library Preparation, Step 02)

Adapter: dilute according to **Table 2** (Page 02).

Cleanup: use 0.6x beads, elute DNA with 22.5 µl of elution buffer, transfer 20 µl supernatant for next step.

#### Step 3: Library Amplification (Refer to 08/ Standard Protocol of Library Preparation, Step 03)

Number of cycles: 12 - 17 is recommended. Adjust according to requirement of library yield.

Cleanup: whether to perform a size-selection depends on the sample situation and data analysis requirements.

▲ No size-selection: purify with 0.9x beads, elute DNA with 22.5 µl of elution buffer. Transfer 20 µl of supernatant to a new Nuclease-free PCR tube and store at -20°C.


▲ Size-selection: size-select with 0.73x / 0.25x beads, elute DNA with 22.5 µl of elution buffer. Transfer 20 µl of supernatant to a Nuclease-free PCR tube and store at -20°C.

## Appendix 2. cfDNA Library Preparation


### Step 4: Library Quality Control

**Library concentration determination:** it is recommended to use fluorescent dye (Equalbit or Picogreen) or QPCR for library quantification.


**Library size distribution detection:** detect by an Agilent 2100 Bioanalyzer.



cfDNA extracted with beads



cfDNA Library (without size-selection)



cfDNA Library (after size-selection)

## Appendix 3. Target Capture Library Preparation

Taking NimbleGen SeqCap EZ capture library workflow as an example, VAHTS Universal Pro DNA Library Prep Kit for Illumina is suitable for the preparation of target capture library.

### Notes

- ◇ "Input DNA" is the DNA added to DNA Damage Repair & End Preparation with volume  $\leq 48 \mu\text{l}$ .
- ◇ The length of input DNA should be between 180 bp - 220 bp. Refer to the instruction of Covaris or other fragmentation instruments for fragmentation parameter.
- ◇ To ensure the quality of library, it is recommended to detect the size distribution (by an Agilent 2100 Bioanalyzer) and concentration (by Equalbit) of input DNA.

### Library Preparation Procedures

#### Step 1: DNA Damage Repair & End Preparation (Refer to 08/ Standard Protocol of Library Preparation, Step 01)

Input DNA: according the sample type, please refer to **Table 1**(Page 02).

#### Step 2: Adapter Ligation (Refer to 08/ Standard Protocol of Library Preparation, Step 02)

Adapter: dilute according to **Table 2** (Page 02). VAHTS DNA Adapters for Illumina (**Vazyme, # N801/N802**) are perfectly compatible with Nimblegen Seqcap EZ. When using capture reagents adapters from other sources, choose adapters according to the Blocking reagents.

Cleanup: by 0.6x beads, elute DNA with 105  $\mu\text{l}$  of elution buffer, transfer 100  $\mu\text{l}$  of supernatant for size-selection (with 0.68x / 0.20x beads), then elute DNA with 22.5  $\mu\text{l}$  of elution buffer and transfer 20  $\mu\text{l}$  of supernatant for next step.

#### Step 3: Library Amplification (Refer to 08/ Standard Protocol of Library Preparation, Step 03)

Number of cycles: according to **Table 3** (Page 04). It is recommended to use the upper limit of the cycle numbers, which is enough for a library yield  $\geq 1 \mu\text{g}$ . With sample pooling before target capture, please ensure that all library yields are  $\geq 1 \mu\text{g}/n$  ( $n$  is the sample number). Under this situation, reduce amplification cycles to increase the library complexity and decrease duplication rates.

Cleanup: purify by 0.9x beads, elute DNA with 22.5  $\mu\text{l}$  of ddH<sub>2</sub>O (instead of Dilution Buffer) and transfer 20  $\mu\text{l}$  of supernatant to a new Nuclease-free EP tube.

#### Step 4: Library Quality Control

Refer to **Nimblegen Seqcap EZ Library SR Users Guide V5.1 (Chapter 4, Step 5)** (Roche document number 06588786001, 09/15)

#### Step 5: Target Enrichment

Refer to **Nimblegen Seqcap EZ Library SR Users' Guide V5.1 (Chapter 5-8)**.



ISO 9001: 2015